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CE 0044

85GADRSOKT - 0919

and extensions should only be used by surgeons, dentists, physicians and dental technicians who have undergone intensive training in implantology. The dentist alone decides on the practical application and the respective design of the implant system as well as on the associated prosthetic equipment and therefore on the detailed use of the product. The dentist will make the appropriate decision depending on the individual situation (indication) and take full responsibility for the chosen approach. The use of inappropriate methods can result in the loss of the implants and in considerable damage to the bone surrounding the product.
Nota: Dental Ratio Systems GmbH offers information and training sessions for physicians, dentists and dental technicians. During intraoral application attention must be paid to the fact that the products are protected against aspiration or falling on the floor.

POTENTIAL SIDE EFFECTS

Some of the potential side effects include: chronic pain, swelling, infection, poor healing of the surgical site, sensitivity of the surrounding tissue, allergic reactions, bleeding, exfoliation, local or systemic infections, uncontrolled bone loss, nerve damage, and temporary or long-term tissue paresthesia. Unsuccessful osseointegration is indicated by infection mobility or bone loss.

Safety, liability and warranty

No structural or surface modifications may be performed on the implant. The use of components not belonging to the system as well as any type of modification may impair the function of the implant system (problems of fit, reduced duration as well as damage of the implant up to fracture) and precludes any guarantee or replacement by the manufacturer or DRS. This applies in particular to other application procedures that have not been recommended. As the processing and use of DRS products is outside our control, the user shall be solely responsible. We do not accept liability for any damage caused due to this.
Advice regarding the use of our products is given verbally, in writing, via electronic media or through demonstrations. Our information is based on the state of science and technology as it is known to us at the time of product distribution. This does not release the user from the duty to personally check the product for suitability with regard to the respective indication and to take part in further training covering the field of dental-alevolar surgery/implantology. It is essential to observe the color coding and/or labels. Thermal damage caused by rotating tools in the jaw must be avoided. (Use training; working with low speeds and sufficient cooling).

Notes

After implantation, the implants should be surrounded on all sides by at least 1.0 mm of residual bone. If this is not the case (thickness of the bone wall is less than 1.0 mm) or the bone lamella is non-existent), a bone augmentation procedure is indicated. Restorative measures must not be introduced until the OSSEOINTEGRATION of the installed implant has succeeded.
DRS implants can be used for all indications of oral, endosteal implants in the upper and lower jaws, for functional and aesthetic oral rehabilitation of fully and partially edentulous patients. The success of the implant is largely dependent on the prosthetics. Intensive communication between dentist and dental technician, careful pre-prosthetic planning and the involvement of the patient are important requirements for the success of the implant. Prosthetic applications can include individual crowns, bridges, partial or total prostheses, which are connected with the DRS implants via the respective elements.

Specific instructions for small diameter implants

Due to the reduced mechanical stability, Ø 3.3 mm implants are used only in cases with low mechanical load.
Warning: Small diameter implants and angled abutments are not recommended for the posterior region.

Special indications for Tissue Level implants with reduced diameter

- Jaws with teeth: 4 implants with a bar for primary connection
- Jaws with some teeth: in the case of implant-borne fixed restorations that are combined with Ø 4.1 mm or Ø 4.8 mm implants and whose superstructure has primary splinting.

Special contraindications for Tissue Level implants with reduced diameter

- Single canine tooth restorations and single tooth restorations in the posterior area
- Use with retentive anchors, magnets or telescope crowns, alligator abutments

Specific instructions for implants with a length of 6 mm

(Warning: 6 mm length implants are not approved for use in Canada and USA)
Because of the reduced anchorage surface area in the bone, the implants with a length of 6 mm are suitable for use with the following indications:

- Edentulous jaw: as an auxiliary implant for implant-supported fixed bar constructions to support full dentures.
- Jaws with some teeth: as an auxiliary implant combined with longer implants whose superstructure has splinting.

Specific instructions for Wide Platform Implants

These implants are indicated for dental replacement in the molar region of the mandible or maxilla where the width of the alveolar ridge is at least 6.8 mm. Clearance of at least 5.0 mm from the center of the implant to the next natural tooth and a medio-distal distance of at least 9.5 mm at bone level are required. Due to the broader axial profile of the Wide Platform implants, reconstructions with superior aesthetics and stability can be made in the molar region than on Regular Platform implants.

Specific instructions for Bone Level Tapered Design Implants

These implants can be used for immediate or early implantation after loss or extraction of natural teeth. The endosseous implant body has a cylindrical shape in the coronal region and is conical in the apical region. The coronal cylindrical part of the implant is relevant for the stated intraosseous diameters.

Contraindications

When selecting implant patients, it is important to observe the general medical as well as local, absolute and relative contraindications for dental surgery. Contraindications include serious internal diseases, bleeding disorders, metabolic disorders, bone-related metabolic disorders (i.e. osteoporosis), reduced immune defense, blood clotting disorders and allergies to the used materials (titanium Grade 4). Relative contraindications are bruxism, psychiatric conditions, use of medication, consumption of drugs, alcohol and tobacco, as well as insufficient bone quantity and defective bone quality and also poor oral hygiene and unhealthy remaining teeth. In the case of insufficient bone quantity, we recommend the use of currently available augmentation procedures and corresponding literature. An insufficient number, the inadequate length or the faulty positioning of implants can lead to biomechanical strain resulting in loss of implants, fatigue fractures of implants, abutments or prosthetic screws. This should be considered especially for prosthetic restorations with retentive anchors, an adverse force distribution may lead to a breakage of prosthetic parts.

Treatment planning

The surgical plan of implants needs to be based on comprehensive patient evaluation, careful pre-surgical diagnostics and documented therapy planning. Mistakes in patient evaluation, pre-surgical diagnostics and therapy planning can lead to a premature loss of the implant.

Healing phase

Manual osseointegration tends to take three months, maxillary healing generally takes six months. Depending on bone quality and anatomy, the lengths of these periods can vary. After the healing phase and gingival shaping (healing abutment), the prosthesis can be applied.

Implant installation procedure

It is also essential to treat the hard soft-tissue with due care in order to ensure optimal conditions for the osseointegration period of the osteoimplant system. The intervention should be carried out as atraumatically as possible and requires a high degree of precision care. The parallel drilling should be carried out using suitable sharp-edged drills and the maximum rotation speed indicated below. The operator must ensure that the tip of the drill is sufficiently cooled to avoid thermal trauma, which can prevent successful osseointegration of the implant. We recommend the use of physiological sterile sodium chloride solution as cooling fluid. Insertion of the implant as well as the rough cutting of a thread should occur at a maximum of 35 Nm (regardless of whether this is done mechanically or manually).

Warning: The primary stability of the implant immediately after installation must be ensured. The intraosseous implant surface must be completely covered by bone or bone substitute during surgery. The final alignment of the implant must be checked after the prosthetic restoration. It is recommended to align the markings on the transfer part (blue) or- facially. Instruments and tools must be secured to avoid unintentional aspiration by the patient.

Notes on specific contact with OKTAGON® implants

- Avoid any contact with or contamination of the implant.
- Please note that the apical part of the drill tips of the spiral drills are overlong, i.e. the preparation depth of the implant bed does not correspond to the insertion depth of the implant.
- This must be taken into account when selecting the implant length.

- Patient evaluation/preoperative treatment planning: It is essential that these steps precede the surgical phase as insufficient planning may lead to a loss of the implant.
- Incision: The gingiva is opened using a scalpel. This is followed by creating a mucoperiosteal flap to reveal the osseous structures.
- Implant positioning bur (max. 800 min^{−1}): The implant site is determined according to the location defined during previous planning.
- Ø 2 mm twist drill (max. 800 min^{−1}): Pre-drilling is performed at the punched site taking into account and controlling the axis of the implant. The planned insertion depth can be determined via the depth markings on the drill. This applies to all twist drill holes. Placing of the parallel pin serves for checking the axis of the implant and, if applicable, parallelism. In addition, the depth can be checked using the depth measuring gauge which is available in different diameters.
- Twist drill: For the Conical implant and the Bone Level Tapered implant, system-specific twist drills are used to prepare implant bed.
Important: After pre-drilling with the Ø 2.2 mm twist drill for Conical implants, use only one system-specific twist drill corresponding to the implant length.
 - Ø 2.8 mm twist drill (max. 600 min^{−1}): Enlargement of the implant bed to Ø 2.8 mm with check of the implant axis, minor directional corrections are still possible here. This is followed by thread-cutting for the OKTAGON® implant Ø 3 mm.
 - Ø 3.25 mm twist drill (max. 500 min^{−1}): Enlargement of implant bed to Ø 3.25 mm for insertion of the implant. This is followed by thread-cutting for the OKTAGON® implant Ø 3.75 mm.
 - Ø 3.5 mm twist drill (max. 500 min^{−1}): Enlargement of implant bed to Ø 3.5 mm. In this case, engraving with Ø 3.25 mm twist drill is not mandatory. This is followed by thread-cutting for the OKTAGON® implant Ø 4.1 mm.
 - Ø 4.2 mm twist drill (max. 400 min^{−1}): Enlargement of implant bed to Ø 4.2 mm. This is followed by thread-cutting for the OKTAGON® implant Ø 4.8 mm.
- Profile drill (max. 300 min^{−1}): The profile drill is used for fire preparation of the implant bed and prepares this specifically for the different types of implants. The implant types Tapered Design, Wide Platform, Bone Level and Conical require implant bed preparation with the profile drill.
- Thread-up (max. 15 min^{−1}): The OKTAGON® is a non-selecting implant system. In the case of compact bone, as (implant-specific) thread needs to be cut as a rule. For Ø 4.8mm implants, thread-up is recommended for all bone classes.

Notes on removing implant and disconnecting transfer component (blue) with a screw:

The sterile packaging is only to be opened in the treatment room directly prior to implantation. The implant is supplied in a vial with a titanium seal, whereby the label with data on the OKTAGON® implant and the LOT No. is transferred to the patient file. The adapter is attached to the screw of the blue transfer component. The implant is removed from the plastic vial by turning to the left. Attention: Please hold vial horizontally when removing the implant!

Prior to using the insertion and check the correct functioning and retaining force of the instrument and ensure that the implant is gripped securely.

- Insertion: After removing the OKTAGON® implant from the vial with the help of the inserion aid, the implant is inserted (manually or automatically) carefully into the prepared implant bed. A torque of 35Nm should not be exceeded during insertion as this could damage the bone, thus increasing the risk of implant loss. To reduce this risk, the thread cutter should be used beforehand. After the implant has been inserted, it is secured against twisting on the interface of the blue transfer component with a locking key, and the transfer component disconnected from the implant by turning the adapter to the left.

Wound closure: The OKTAGON® implant system can be used for both transgingival (open) or subgingival (submerged) healing. In the case of submerged healing, the inner thread of the implant is closed with the aid of a cover screw, a suture procedure is indicated. Numerous items are available in the product range of transgingival healing, which form part of the soft tissue during open healing. To this purpose the wound margins are sutured atraumatically tight around the restoration. The sutures must not impair blood supply to the vessels of the gingival margin.

Re-opening surgery (subgingival healing): Following the healing phase, the gingiva is opened surgically at the implantation site, the cover screw is removed and replaced with a suitable abutment.

15. Prosthetic restoration: Using the impression coping now allows impression-taking and a model to be produced using laboratory analogues. The restorative prosthesis can be fabricated and screwed to the implant. A maximum torque of 35Nm is recommended for the screw connection. For detailed information please contact DRS International GmbH.

General information / MRI compatibility

The patient must receive sufficient, detailed instructions about the oral hygiene procedures needed. It is recommended to include the patient in a regular recall system. The dental Implant OKTAGON® has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Dental Implant OKTAGON® in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Packaging and sterility

OKTAGON® implants sold by the company Dental Ratio Systems GmbH are supplied in sterile form (gamma-ray sterilization, minimum strength 25kGy). They are sold in special plastic tubes, fastened on a titanium support with a pre-assembled insertion system. The tube with the implant is double-wrapped in special blister packaging. The label can be peeled off and used to enter the implant data in the patient file (lot no.). OKTAGON® implants are designed to be single use only! The implant is sterile as long as the packaging is unopened and undamaged. If the packaging is damaged or out of the inner wrapping is broken, there is a risk of contamination, i.e. the implant must not be used or re-sterilized. OKTAGON® implants must not be used beyond the expiry date of guaranteed sterilization (hourglass symbol). Unprotected touching of the implant by the user must be avoided (use of protective gloves is recommended). All parts of the OKTAGON® implant system must be stored in a dry location at room temperature and, if possible, should not be exposed to direct sunlight. The re-use of single-use products is associated with a high risk of infection, and risk-free functional safety (e.g. regarding stability or fitting accuracy) cannot be guaranteed.

Storage, disinfection, cleaning and sterilization of abutments, components and tools

Apart from the implants themselves, all other parts of the OKTAGON® implant system (abutments and instruments) are supplied non-sterilized. Prior to the initial use on the patient and, if possible, immediately after each use, all components and instruments of the OKTAGON® implant system that have been supplied non-sterilized, must be disinfected, cleaned and sterilized – as far as this is necessary and possible regarding their intended use. If the sterilization method is not specified, then use the disinfection method used by the manufacturer. Disinfectants used for this purpose are alcohol, peracetic acid, glutaraldehyde, hydrogen peroxide, or steam. Only use disinfectants if the products should be stored in their original packaging at room temperature and protected from dust and humidity. Failure to clean and sterilize the components according to the instructions may result in harmful infections. All suitable components can be sterilized in accordance with DIN EN ISO 17665 - Sterilization of health care products - Moist heat and STPS (30 min at 121 °C (250 °F) gravity - displacement with 20 mins drying time and/or 4 minutes for 132° (270° F) / 134° (273° F) fractionated vacuum with 20 minutes drying time). The STPS need to be wrapped for sterilization, use sterilization pouches for wrapping. Please observe the manufacturer's instructions for sterilization. Prior to sterilization, components must be adequately disinfected and cleaned. Disinfection can be carried out manually or mechanically. When instruments are very dirty, use ultrasonic cleaning. After cleaning, the instruments must be rinsed thoroughly with water and dried immediately. The instruments must not come into direct contact with other instruments during the cleaning/disinfection process in order to avoid damage. Always wear protective gloves when handling dirty instruments. Sterilization temperatures in excess of 134°C (273° F) must be avoided. Exceeding the maximum storage duration must not be exceeded. All rotating instruments are, in principle, suited to multiple use. The dentist/expert alone is responsible for any further use of an instrument or for the frequency of use, and will base this decision on the respective case as well as on the potential wear and tear of the instrument. When in doubt, the instruments should be refused and replaced in good time. All structural components are, in principle, suited to multiple use. However, the dentist, i.e. the user, alone decides whether, based on the respective case as well as on the potential wear and tear leading to material fatigue of some components, multiple use makes sense. Cutting instruments can be reused a maximum of 10 times (1 use = 1 implant inserted) as long as they are adequately handled and are undamaged and uncontaminated. Any other use beyond this use or the use of damaged or contaminated instruments is prohibited.

RPM recommendations for rotating instruments

The best working results are achieved when the user observes the instrument-specific rotational speed recommendations. Long and pointed instruments exceeding the maximum admissible number of revolutions tend to develop vibrations that can lead to the destruction of the instrument. When the diameter of working parts exceeds the strength of the shaft, an excessive number of revolutions may result in strong centrifugal forces that lead to distortions of the shaft and/or to the instrument breaking. For this reason, the admissible maximum number of revolutions must never be exceeded. You will find information on the recommended rotational speed and the maximum admissible rotational speed in the data provided by the manufacturer. The non-observation of the admissible maximum number of revolutions results in a heightened safety risk.

Please follow general application and safety instructions for MEISINGER products in the medical area and also the advice for re-processing (cleaning, disinfection and sterilization) of medical devices from Hager & Meisinger GmbH.

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Systèmes d'implant OKTAGON® (Tissue Level : Regular Platform, Wide Platform, Tapered Design, Conical/Bone Level : Narrow Connect, Regular Connect, Tapered), implants, piliers, composants/accessoires chirurgicaux, de restauration et de dentisterie

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